1 Einführung

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Patent Documents: General Information Warranty Safety and Handling

1 Introduction

1.1 Copyright and intellectual property
The copyright for these assembly and use instructions is held solely by the manufacturer. Additionally, all rights are reserved, particularly with regard to patent grant or utility model registration. Violations of the above-mentioned provisions will result, where applicable, in the obligation to pay compensation for damages.
Intellectual Property – Patent numbers
U.S. Patent No.: 8529484
CA Patent No.: 2692534
Patent pending in Europe

1.2 Foreword
We would like to thank you for your confidence in TurboMed and the FS3000. This orthosis was designed by TurboMed Orthotics Inc., located in the province of Quebec, in Canada. This new-generation orthotic solution was developed specifically for people with foot drop.
The FS3000 was developed to offer you unrivalled performance and comfort.
We urge you to read this document carefully and address any questions and/or comments to your authorized retailer before using or assembling your new orthosis. In order to ensure that the equipment works properly and safely, you must strictly follow all of the instructions provided in this document and observe the specified usage restrictions. Do not hesitate to contact your authorized retailer.
If you have any questions about the way your device is working or if you have any questions, do not hesitate to contact your authorized retailer.
To contact the customer service department of TurboMed Orthotics in Canada, please visit: www.TurboMedOrthotics.com
You'll find a complete list of distributors by country on our website.

1.3 Use
Intended use: This orthosis is intended SOLELY as an orthotic fitting for the lower limb and has been designed ONLY as a foot lift (droptop) device.
It is important to use it in accordance with the indications.

1.4 Safety symbols used in this guide
The symbols below are used in this guide to identify the warnings and useful advice that you should observe. It is very important to read and thoroughly understand them before using our product for the first time.

Warning: this symbol means that there is a potential risk. Failure to heed a warning may result in bodily harm or injury.

Notice: recommendations or instructions should be taken seriously.

Product-related caution: concerns a potential material risk. Failure to heed a product caution may result in significant material damage.

2 Presentation of the product

Item	Description
1	MAIN POSTERIOR SECTION (Y" SHAPE)
2	FRONT HARNESS (U" SHAPE)
3	BACK HARNESS
4	LACE CLIP
5	CALF SECTION
6	ALUMINUM JUNCTION BLOCK

3 Indications and contraindications

3.1 Indications
The orthosis is intended for people who have:
Weakness in the muscles that lift the foot, presenting no spasticity to mild spasticity (post. back, hip or knee surgery; traumatic injury to the lower limbs; a stroke; traumatic brain injury; multiple sclerosis; Charcot-Marie-Tooth disease; peroneal paralysis; neuromuscular atrophy, etc.).
No significant deformities of the foot that require a visit to a physician to verify whether the orthosis is appropriate for the client or not.
Signs of this product must always be determined by the doctor and/or prescribing physician.

3.2 Contraindications
The orthosis was not designed for cases of:
-moderate to severe spasticity of the lower limb;
-leg ulcers in the posterior part of the calf;
-severe balance problems;
-moderate to severe foot deformities;
-certain pathologies that require a visit to a physician to verify whether the orthosis is appropriate for the client or not.

4 The product

4.1 Usage restrictions
Use of the orthosis is not recommended to prevent contractions (shortening of muscles), etc.
Avoid going near a source of high heat, such as a fire, a heating element, embers, flames, etc.
Avoid using an orthosis that belonged to someone else or an orthosis that you believe has been used previously.

4.2 Warning and general safety instructions (client)
Information:
-please read this document in its entirety before using the product;
-heed the safety instructions in order to avoid injury or damage to the product;
-keep this document in a safe place.

4.3 Warning and general safety instructions (professional)
Information:
-please read this document in its entirety before using the product;
-need the safety instructions in order to avoid injury or damage to the product;
-show users how to properly use your product and inform them accurately of the safety instructions;
-keep this document in a safe place.

4.4 Caution
⚠
This product is for single patient use only.
Stop using the product regularly if there are changes in fit, stiffness or weight.
Check the product regularly to see if it still works properly.
Stop using the product if there are changes in, or loss of, functionality and have it inspected by your authorized retailer.

Observe the legal and insurance regulations relating to operating motor vehicles of any type and have your driving licence tested by the competent authority.
In general, the physician will determine the daily duration of use and the total period of use.
Only qualified employees are authorized to perform the initial fitting and the first application of the product.
Do not use the orthosis until you have received the necessary training.
Inform the patient that it is important to see a physician immediately if unusual changes to the body are observed (ex., the development of a wound, worsening of pain, etc.).
Make sure that the orthosis is applied correctly and correctly.
Inform the patient to check, prior to each use, that the product works properly and that it shows no signs of wear or damage.
Direct the patient's attention to the fact that he/she CANNOT continue to use the product if it or one of its components is worn or damaged (ex., cracks, missing part, warping, inappropriate shape of the orthosis, etc.).

5 Application and removal of the product
5.1 Application
1 Attach the lace clip to an appropriate shoe.
2 Open the calf band.
3 Install the product onto the shoe.
4 Put on the shoe and orthosis, pushing back on the orthosis.
5 Close and tie the shoe.
6 Lose the calf band.

5.2 Removal
1 Open the calf band.
2 Untie your shoe.
3 Remove your foot from the shoe.

6 Care and cleaning
Clean the rigid parts of the product only with mild soap and water or a damp cloth. The so-called "soft" parts, such as the calf pad, the band, etc., should be cleaned by hand using mild soap, then air-dried. Avoid any source of direct heat (ex., solar radiation, a radiator, fire, etc.).

7 Garanties

7.1 Garantie und assistance
Die FS3000 orthosis was manufactured by TurboMed Orthotics Inc. TurboMed guarantees its orthosis for a period of two years (2 years) from the date of purchase. If a product proves to have material or manufacturing defects, TurboMed will repair it or replace it at its own expense if sufficient proof is provided—namely, a receipt—and if there is proof that the date of the service request is within the guarantee period. The use of unauthorized replacement parts will cancel this guarantee. Contact your authorized retailer or TurboMed for additional information in this regard.

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7.2 Spezialanfertigungsgarantie
We are so confident in the quality of our product that we offer to refund the cost of your orthosis if you are not entirely satisfied (see above) within 2 months of your purchase.
turboMed, its distributors and their authorized retail agents undertake jointly and severally to refund in full the amount paid to the client if he/she is not satisfied with the product within two months (2 months) of the date of purchase of the TurboMed orthosis. The client and/or authorized retailers and/or distributors will be solely responsible for the shipping expenses related to sending the orthosis back to us. However, prior authorization must be requested.

8 Legal information
All legal requirements are subject to the national legislation of the country in question and therefore may vary depending on where the product is used.

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7.1 保修和协助
Die FS3000 orthosis was manufactured by TurboMed Orthotics Inc. TurboMed guarantees its orthosis for a period of two years (2 years) from the date of purchase. If a product proves to have material or manufacturing defects, TurboMed will repair it or replace it at its own expense if sufficient proof is provided—namely, a receipt—and if there is proof that the date of the service request is within the guarantee period. The use of unauthorized replacement parts will cancel this guarantee. Contact your authorized retailer or TurboMed for additional information in this regard.

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